

Cellix Ltd.

Inish Device

Portable Hematology Analyzer



Inish Device Advantages

- Fast, reliable, label-free Cell and Particle analysis:
- → detection, characterization & viability
- Portable, low cost device; works off 12V car battery → new applications possible.
- Technology as sensitive as traditional flow cytometry; benchmarked against gold standards: BD's Accuri and Sysmex XP-300. CVs < 3%.





RUN 1

RUN 2


Inish Device: Haematology Analyzer – Validation Data for RBC Counts

Benchmarking Inish against BD Accuri:

- Flow Rate: 35µl/min
- 10 samples at 1:4,000 serial dilution (1:20 dilution from stock)
- Run simultaneously
- Analysed using Inish software and ACCURI software

Benchmarking of Inish vs ACCURI: RBC counts: thousands per ul

RUN 3



RUN 4


RUN 5




Inish Device: Haematology Analyzer – Validation Data for WBC viability

Benchmarking Inish against BD Accuri:

- Sample:
 - 1. 70% live + 30% dead (ethanol added)
 - 2. 100 % Dead
- Flow Rate: 35µl/min
- Staining:
 - o Calcein (live)
 - Ethidium homodimer 1 (dead)
- Run simultaneously
- Analysed using Inish software and ACCURI software (FL1 & FL3 based thresholding)







For more information, please contact us at:

Contact Details:

Tel: +353 (1) 4500 155

Email: info@cellixltd.com

Web: www.cellixltd.com